ENGLISH

NX3 Universal Resin Cement System

NX3 is a permanent resin cement system that offers innovative chemistry for unmatched esthetics, delivery system versatility and compatibility with both total-etch and self-etch adhesives. Delivery system choices include an automatic dual syringe dual-cure cement (base/ catalyst) and a single-syringe light-cure (base only) cement for multiple venter units where flexible work time is needed. NX3 dual-syringe can be used for all indirect applications including veneers and metal-based restorations, without the need of an activator for the adhesive. A proprietary amine-free initiator system provides excellent color stability for long-term esthetics. One set of try-ins gets match both light-cure and dual-cure resin cements. NX3 contains nano filler technology for improved dispersing, handling, cleanup, and polishing properties.

Indications for use:

Cementation of veneers, inlays, onlays, crowns, bridges, and posts. Adhesive bonding of amalgam restorations. Core-buildup material.

Restorative materials:

Indicated for cementation of ceramic, porcelain, resin, metal-based materials and CAD/CAM blocks.

Note: *When a total-etch technique is desired, the use of OptiBond™ Solo Plus** adhesive with NX3 resin cement is recommended. When a self-etch technique is desired, the use of OptiBond XTR*** is recommended.*

** o OptiBond S

*** o OptiBond Versa

I. Fit Evaluation of Restoration, Try-in, and Restoration Surface Preparation

- Remove temporary restoration and thoroughly clean the tooth (teeth) with fluoride-free prophy paste. Rinse thoroughly with water and lightly dry with air.
- Evaluate fit. Make any necessary adjustments.
- To evaluate shade (for composite and ceramic restorations), place a thin layer of the selected water-soluble try-in paste on the restoration. The try-in paste matches both the single-syringe light-cure cement, and the dual-syringe dual-cure cement. Using light pressure, completely seal the restoration. Remove any excess try-in paste.

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY OR FITNESS FOR A PARTICULAR PURPOSE.

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- Work time is based on material at 23°C. Clean-up time is based on material at 37°C. This time may vary based on storage conditions, temperature, humidity, and age of product.
For the dual-cure cement, bleed cartridge before initial use.
- Both single-syringe light-cure resin cement and dual-syringe, dual-cure resin cement (ideal for 1-2 units) can be used for cementing the veneers. For thicker veneers, a dual-cure resin cement is recommended.
- Dispense cement directly into the internal surface of the veneer. Place the veneer gently on the tooth, allowing the cement to slowly flow from all margins. Spot cure the veneer into place on the facial surface away from the margins using a small diameter light guide for 10 seconds. After removing the excess cement, light cure all surfaces for a minimum of 20 seconds* per surface.
- B. Crowns, bridges, inlays, onlays, and metal-based restorations**
Note: On the prep, avoid pooling of adhesive before light curing of adhesive. Excess adhesive can be removed with a dry applicator brush before light-curing of the adhesive.
Apply the dual-cure cement to the restoration or the prep. Seal the restoration gently onto the preparation allowing the cement to flow from all sides. Remove excess cement. Light cure all surfaces for a minimum of 20 seconds* per surface.
† Excess cement cleanup – Excess cement is best removed in its gel state with a scaler or explorer. Gel state can be achieved by tack curing with a light for 1-2 seconds or allowing the cement to self-cure for 2-3 minutes after application.

Post Bonding and Core buildup
Note: On the prep, avoid pooling of adhesive before light curing of adhesive. Excess adhesive can be removed with a dry applicator brush or absorbent paper point before light-curing of the adhesive.
1. Prepare post space. Size and fit post.
2. Place adhesive on the post in a thin even layer, lightly air thin if necessary. Light-curing of the adhesive is optional.
3. Apply the dual-cure cement onto the post and/or into the post preparation, seal the post, and vibrate the post slightly to allow the composite of trapped air.
4. When you are satisfied that the post is properly sealed, remove all excess cement. Light cure all surfaces for a minimum of 20 seconds* per surface.
5. Proceed with core buildup and/or preparation according to manufacturer's instructions for use.

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